

## Rabbit Anti-Cat IgG (H&amp;L) - AF568

Catalog No: #L35324



Package Size: #L35324-1 100 ul #L35324-2 500 ul #L35324-3 1 ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rabbit Anti-Cat IgG (H&L) - AF568
Brief Description	Conjugated Secondary Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, IH, IF, FC
Species Reactivity	Ct
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Cat IgG. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Cat IgG
Conjugates	AF568
Formulation	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide

## Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only