

Mouse Anti-Human Ig lambda chain - AF568

Catalog No: #L35352



Package Size: #L35352-1 100 ul #L35352-2 500 ul #L35352-3 1 ml

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Anti-Human Ig lambda chain - AF568
Brief Description	Conjugated Secondary Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, IH, IF, FC
Species Reactivity	H
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Human Ig lambda chain. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Human Ig lambda chain
Conjugates	AF568
Formulation	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide

Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only