

## Rabbit Anti-Guinea Pig IgG (H&amp;L) - AF633

Catalog No: #L35416



Package Size: #L35416-1 100 ul #L35416-2 500 ul #L35416-3 1 ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rabbit Anti-Guinea Pig IgG (H&L) - AF633
Brief Description	Conjugated Secondary Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, IH, IF, FC
Species Reactivity	Gp
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Guinea Pig IgG. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Guinea Pig IgG
Conjugates	AF633
Formulation	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide

## Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only