

## Mouse Anti-Duck IgY (H&amp;L) - AF633

Catalog No: #L35429



Package Size: #L35429-1 100 ul #L35429-2 500 ul #L35429-3 1 ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Anti-Duck IgY (H&L) - AF633
Brief Description	Conjugated Secondary Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, IH, IF, FC
Species Reactivity	Dk
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Duck IgY. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Duck IgY
Conjugates	AF633
Formulation	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide

## Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only