

Goat Anti-Human IgG F(ab')<sub>2</sub> - AF647

Catalog No: #L35482



Package Size: #L35482-1 100 ul #L35482-2 500 ul #L35482-3 1 ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |   |
|-----------------------|---|
| Product Name          | Goat Anti-Human IgG F(ab') <sub>2</sub> - AF647   |
| Brief Description     | Conjugated Secondary Antibody   |
| Host Species          | Goat  |
| Clonality             | Polyclonal  |
| Purification          | The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.  |
| Applications          | E, IH, IF, FC   |
| Species Reactivity    | H   |
| Specificity           | By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgG F(ab') <sub>2</sub> . No antibody was detected against non immunoglobulin serum proteins. |
| Immunogen Description | Human IgG F(ab') <sub>2</sub>   |
| Conjugates            | AF647   |
| Formulation           | 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide   |

## Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only