

## Goat Anti-Human IgG Fab - AF680

Catalog No: #L35529



Package Size: #L35529-1 100 ul #L35529-2 500 ul #L35529-3 1 ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |  |
|-----------------------|--|
| Product Name          | Goat Anti-Human IgG Fab - AF680  |
| Brief Description     | Conjugated Secondary Antibody  |
| Host Species          | Goat   |
| Clonality             | Polyclonal   |
| Purification          | The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.  |
| Applications          | E, IH, IF, FC  |
| Species Reactivity    | H  |
| Specificity           | By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgG Fab. No antibody was detected against non immunoglobulin serum proteins. |
| Immunogen Description | Human IgG Fab  |
| Conjugates            | AF680  |
| Formulation           | 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide  |

## Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only