

## Mouse Anti-Monkey IgG (H&amp;L) - AF680

Catalog No: #L35564



Package Size: #L35564-1 100 ul #L35564-2 500 ul #L35564-3 1 ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |   |
|-----------------------|---|
| Product Name          | Mouse Anti-Monkey IgG (H&L) - AF680   |
| Brief Description     | Conjugated Secondary Antibody   |
| Host Species          | Mouse   |
| Clonality             | Monoclonal  |
| Purification          | The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.  |
| Applications          | E, IH, IF, FC   |
| Species Reactivity    | Mk  |
| Specificity           | By immunoelectrophoresis and ELISA this antibody reacts specifically with Monkey IgG. Do not cross reactivity with human IgG. No antibody was detected against non immunoglobulin serum proteins. |
| Immunogen Description | Monkey IgG  |
| Conjugates            | AF680   |
| Formulation           | 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide   |

## Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only