

Chicken Superoxide Dismutase [Cu-Zn] (SOD1) ELISA Kit



Catalog No: #EK20443

Orders: order@signalwayantibody.com

Package Size: #EK20443-1 48T #EK20443-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Chicken Superoxide Dismutase [Cu-Zn] (SOD1) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Chicken (Gallus)
Storage	12 months

Application Details

Detect Range:15.625-1000ng/mL

Sensitivity: 6.52ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 µL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only