

# Chicken NADH Dehydrogenase, Quinone 1 (NQO1) ELISA Kit



Catalog No: #EK20403

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20403-1 48T #EK20403-2 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Chicken NADH Dehydrogenase, Quinone 1 (NQO1) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Chicken (Gallus)
Storage	12 months

## Application Details

Detect Range:0.31-20ng/mL

Sensitivity: 0.14 ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 µL

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only