

# Chicken Reactive Oxygen Species(ROS) ELISA Kit

Catalog No: #EK20447



Package Size: #EK20447-1 48T #EK20447-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Chicken Reactive Oxygen Species(ROS) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Chicken (Gallus)
Storage	12 months

## Application Details

Detect Range:31.25-2000pg/mL

Sensitivity:12.8pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 µL

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only