

NKX3-2 Antibody

Catalog No: #46626

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NKX3-2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Human Mouse
Specificity	The antibody detects endogenous levels of total NKX3-2 protein.
Immunogen Description	Synthetic peptide of human NKX3-2
Target Name	NKX3-2
Other Names	SMMD; BAPX1; NKX3B; NKX3.2
Accession No.	Swiss-Prot:P78367 NCBI Gene ID:579NCBI Protein:NP_001180
Uniprot	P78367
GeneID	579;
Concentration	0.5 mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1: 25-100

Images



The image on the left is immunohistochemistry of paraffin-embedded Human gastric cancer tissue using 46626(NKX3-2 Antibody) at dilution 1/25



The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using 46626(NKX3-2 Antibody) at dilution 1/25

Background

This gene encodes a member of the NK family of homeobox-containing proteins. The encoded protein may play a role in skeletal development.

Note: This product is for in vitro research use only