

NLGN1 Antibody

Catalog No: #46627

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | NLGN1 Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | Antigen affinity purification |
| Applications | IHC |
| Species Reactivity | Hu Rt |
| Specificity | The antibody detects endogenous levels of total NLGN1 protein. |
| Immunogen Type | peptide |
| Immunogen Description | Synthetic peptide corresponding to internal residues of human NLGN1 |
| Target Name | NLGN1 |
| Other Names | NL1 |
| Accession No. | Swiss-Prot:Q8N2Q7 NCBI Gene ID:22871NCBI Protein:NP_055747 |
| Uniprot | Q8N2Q7 |
| GeneID | 22871; |
| Concentration | 2mg/ml |
| Formulation | Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol. |
| Storage | Store at -20°C |

Application Details

Immunohistochemistry: 1: 40-200

Images



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using 46627(NLGN1 Antibody) at dilution 1/50, on the right is treated with synthetic peptide. (Original magnification: x200)

Background

This gene encodes a member of a family of neuronal cell surface proteins. Members of this family may act as splice site-specific ligands for beta-neurexins and may be involved in the formation and remodeling of central nervous system synapses.

Note: This product is for in vitro research use only