

NOL8 Antibody

Catalog No: #46629

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NOL8 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total NOL8 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide corresponding to internal residues of human NOL8
Target Name	NOL8
Other Names	NOP132; C9orf34; bA62C3.3; bA62C3.4
Accession No.	Swiss-Prot:Q76FK4 NCBI Gene ID:55035NCBI Protein:NP_001243323
Uniprot	Q76FK4
GeneID	55035;
Concentration	2mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1: 30-150

Images



The image on the left is immunohistochemistry of paraffin-embedded Human gastric cancer tissue using 46629(NOL8 Antibody) at dilution 1/40, on the right is treated with synthetic peptide. (Original magnification: x200)

Background

NOL8 binds Ras-related GTP-binding proteins (see MIM 608267) and plays a role in cell growth (Sekiguchi et al., 2004 [PubMed 14660641]).

Note: This product is for in vitro research use only