

## RCOR3 Antibody

Catalog No: #46656

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	RCOR3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total RCOR3 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic protein corresponding to residues near the C terminal of human RCOR3
Target Name	RCOR3
Accession No.	Swiss-Prot:Q9P2K3NCBI Gene ID:55758NCBI Protein:BC031608
Uniprot	Q9P2K3
GeneID	55758;
Concentration	0.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C

## Application Details

Immunohistochemistry: 1: 20-100

## Images



The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using 46656(RCOR3 Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: x200)

## Background

May act as a component of a corepressor complex that represses transcription.

---

Note: This product is for in vitro research use only