

# Human cutaneous lymphocyte-associated antigen,CLA ELISA Kit

Catalog No: #EK0192

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK0192 1x96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human cutaneous lymphocyte-associated antigen,CLA ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Hu
Target Name	chemokine (C-C motif) ligand 27
Other Names	ALP, CTACK, CTAK, ESKINE, ILC, PESKY, SCYA27, CC chemokine ILC IL-11 Ralpha-locus chemokine cutaneous T-cell attracting chemokine skinkine small inducible cytokine A27 small inducible cytokine subfa,CCL27
Accession No.	Swiss-Prot: Q9Y4X3
Uniprot	Q9Y4X3
GeneID	10850;

## Application Details

Detect Range: 0.313ng/ml-20ng/ml

Sensitivity: 0.043ng/ml

Sample Type:

Sample Volume: 50-100ul

Assay Time: 1-5h

Detection wavelength: 450 nm

Note: This product is for in vitro research use only