

Nestin Polyclonal Conjugated Antibody

Catalog No: #C29097



Package Size: #C29097-AF350 100ul #C29097-AF405 100ul #C29097-AF488 100ul
 #C29097-AF555 100ul #C29097-AF594 100ul #C29097-AF647 100ul
 #C29097-AF680 100ul #C29097-AF750 100ul #C29097-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	Nestin Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human Nestin (NP_006608.1).
Other Names	NES;Nbla00170;nestin;Nestin
Accession No.	GeneID:10763Swiss Prot:P48681
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	177kDa
SDS-PAGE MW	180kDa
Concentration	1.0mg/ml
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.157.

Application Details

WB □ 1:500 - 1:1000 IHC □ 1:50 - 1:100 IF □ 1:50 - 1:100

Background

This gene encodes a member of the intermediate filament protein family and is expressed primarily in nerve cells.

Note: This product is for in vitro research use only