

α -SMA Monoclonal Conjugated Antibody

Catalog No: #C40482



Package Size: #C40482-AF350 100ul #C40482-AF405 100ul #C40482-AF488 100ul
 #C40482-AF555 100ul #C40482-AF594 100ul #C40482-AF647 100ul
 #C40482-AF680 100ul #C40482-AF750 100ul #C40482-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	α -SMA Monoclonal Conjugated Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous α -SMA protein.
Immunogen Type	peptide
Immunogen Description	Synthetic Peptide
Target Name	α -SMA
Accession No.	Swiss-Prot: P62736NCBI Gene ID: 59
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
SDS-PAGE MW	42kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:5000-50000

Immunohistochemistry: 1:1000-2000

Note: This product is for in vitro research use only