

COL1A2 Polyclonal Conjugated Antibody

Catalog No: #C40771



Package Size: #C40771-AF350 100ul #C40771-AF405 100ul #C40771-AF488 100ul
 #C40771-AF555 100ul #C40771-AF594 100ul #C40771-AF647 100ul
 #C40771-AF680 100ul #C40771-AF750 100ul #C40771-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	COL1A2 Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Species Reactivity	Hu Ms Rt
Specificity	COL1A2 Polyclonal Antibody detects endogenous levels of COL1A2 protein.
Immunogen Description	Synthesized peptide derived from the N-terminal region of human COL1A2.
Target Name	COL1A2
Other Names	COL1A2; Collagen alpha-2(I) chain; Alpha-2 type I collagen
Accession No.	Swiss-Prot: P08123NCBI Gene ID: 1278
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
SDS-PAGE MW	200kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at 4°C in dark for 6 months

Application Details

Western Blot: 1:500-2000

Immunohistochemistry: 1/100 - 1/300

Immunofluorescence: 1/200 - 1/1000

ELISA: 1/10000

Note: This product is for in vitro research use only