

## Desmin Polyclonal Conjugated Antibody

Catalog No: #C40839



Package Size: #C40839-AF350 100ul #C40839-AF405 100ul #C40839-AF488 100ul  
 #C40839-AF555 100ul #C40839-AF594 100ul #C40839-AF647 100ul  
 #C40839-AF680 100ul #C40839-AF750 100ul #C40839-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Desmin Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Species Reactivity	Hu Ms Rt
Specificity	Desmin Polyclonal Antibody detects endogenous levels of Desmin protein.
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Desmin.
Target Name	Desmin
Other Names	DES; Desmin
Accession No.	Swiss-Prot: P17661NCBI Gene ID: 1674
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
SDS-PAGE MW	55kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at 4°C in dark for 6 months

## Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/10000.

Not yet tested in other applications.

Note: This product is for in vitro research use only