

IL-21 Polyclonal Conjugated Antibody

Catalog No: #C41060



Package Size: #C41060-AF350 100ul #C41060-AF405 100ul #C41060-AF488 100ul
 #C41060-AF555 100ul #C41060-AF594 100ul #C41060-AF647 100ul
 #C41060-AF680 100ul #C41060-AF750 100ul #C41060-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-21 Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Species Reactivity	Hu Ms
Specificity	IL-21 Polyclonal Antibody detects endogenous levels of IL-21 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human IL-21.
Target Name	IL-21
Other Names	IL21; Interleukin-21; IL-21; Za11
Accession No.	Swiss-Prot: Q9HBE4NCBI Gene ID: 59067
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
SDS-PAGE MW	18kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Note: This product is for in vitro research use only