

Rat Brevican (BCAN) ELISA Kit

Catalog No: #EK20348



Package Size: #EK20348-1 48T #EK20348-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|-------------------------------|
| Product Name | Rat Brevican (BCAN) ELISA Kit |
| Brief Description | ELISA Kit |
| Applications | ELISA |
| Species Reactivity | Rt |
| Storage | 12 months |

Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.145ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only