

Histone H2B (Phospho-Ser37) Antibody

Catalog No: #SAB294



Package Size: #SAB294-1 50ul #SAB294-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Histone H2B (Phospho-Ser37) Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide. |
| Applications | Custom antibody |
| Species Reactivity | Ms |
| Immunogen Type | Peptide-KLH |
| Target Name | Histone H2B |
| Other Names | CgPICR_009869 |
| Accession No. | uniprot:A0A3L7HVF4 |
| Calculated MW | 60kDa |
| Formulation | Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |

Application Details

Western blotting: 1:500~1:1000

Note: This product is for in vitro research use only