

MFAP3L (Phospho-Tyr287) Antibody

Catalog No: #SAB256P



Package Size: #SAB256P-1 50ul #SAB256P-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | MFAP3L (Phospho-Tyr287) Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide. |
| Applications | Custom antibody |
| Species Reactivity | Hu |
| Immunogen Type | Peptide-KLH |
| Target Name | MFAP3L |
| Other Names | Microfibrillar-associated protein 3-like; Testis development protein NYD-SP9 |
| Accession No. | uniprot:O75121 |
| Calculated MW | 45kDa |
| Formulation | Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |

Application Details

Western blotting: 1:500~1:1000

Background

May participate in the nuclear signaling of EGFR and MAPK1/ERK2. May have a role in metastasis.

Note: This product is for in vitro research use only