

## Arr2 (Phospho-Ser366) Antibody

Catalog No: #SAB345P



Package Size: #SAB345P-1 50ul #SAB345P-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Arr2 (Phospho-Ser366) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	Custom antibody
Species Reactivity	Fruit fly
Immunogen Type	Peptide-KLH
Target Name	Arr2
Other Names	Phosrestin-1; Arrestin-2
Accession No.	uniprot:P19107
Calculated MW	45kDa
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

## Application Details

Western blotting: 1:500~1:1000

## Background

Probably plays an important role in the photoreceptor transduction.

Note: This product is for in vitro research use only