

## PAQR3 (Phospho-Ser8) Antibody

Catalog No: #SAB524P



Package Size: #SAB524P-1 50ul #SAB524P-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PAQR3 (Phospho-Ser8) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	Custom antibody
Species Reactivity	Hu Ms Rt
Immunogen Type	Peptide-KLH
Target Name	PAQR3
Other Names	Progestin and adipoQ receptor family member III,Raf kinase trapping to Golgi,RKTG
Accession No.	uniprot:Q6TCH7
Calculated MW	36kDa
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

## Application Details

Western blotting: 1:500~1:1000

## Background

Functions as a spatial regulator of RAF1 kinase by sequestering it to the Golgi.

Note: This product is for in vitro research use only