

# UCK2 Rabbit Polyclonal Antibody

Catalog No: #29691

Package Size: #29691-1 50ul #29691-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	UCK2 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human UCK2 (NP_036606.2).
Other Names	UCK2;TSA903;UK;UMPK
Accession No.	Swiss Prot:Q9BZX2GeneID:7371
Calculated MW	12kDa/29kDa
SDS-PAGE MW	29kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using UCK2 at 1:1000 dilution.

## Background

This gene encodes a pyrimidine ribonucleoside kinase. The encoded protein (EC 2.7.1.48) catalyzes phosphorylation of uridine and cytidine to uridine monophosphate (UMP) and cytidine monophosphate (CMP), respectively.

---

Note: This product is for in vitro research use only