

Rabbit Glial Fibrillary Acidic Protein (GFAP) ELISA Kit

Catalog No: #EK20131



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rabbit Glial Fibrillary Acidic Protein (GFAP) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rabbit (<i>Oryctolagus cuniculus</i>)
Storage	12 months

Application Details

Detect Range:0.781-50ng/mL

Sensitivity:0.22ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 uL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only