

# Recombinant Mouse Lymphotactin protein(Xcl1) (Active)

Catalog No: #AP73665

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP73665-1 5ug #AP73665-2 20ug #AP73665-3 100ug #AP73665-4 250ug

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

#AP73665-5 500ug #AP73665-6 1mg

## Description

Product Name	Recombinant Mouse Lymphotactin protein(Xcl1) (Active)
Brief Description	Recombinant Protein
Host Species	E.Coli
Uniprot	P47993
GeneID	16963;
Tag Info	NO-tagged
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>

Note: This product is for in vitro research use only