

Recombinant Human Placenta growth factor protein(PGF) (Active)



Catalog No: #AP73561

Orders: order@signalwayantibody.com

Package Size: #AP73561-1 5ug #AP73561-2 25ug #AP73561-3 100ug #AP73561-4 250ug

Support: tech@signalwayantibody.com

#AP73561-5 500ug #AP73561-6 1mg

Description

Product Name	Recombinant Human Placenta growth factor protein(PGF) (Active)
Brief Description	Recombinant Protein
Host Species	E.Coli
Uniprot	P49763
GeneID	5228;
Tag Info	NO-tagged
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>

Note: This product is for in vitro research use only