

Recombinant Human Fibroblast growth factor 1 protein(FGF1) (Active)



Catalog No: #AP73396

Orders: order@signalwayantibody.com

Package Size: #AP73396-1 10ug #AP73396-2 50ug #AP73396-3 100ug #AP73396-4 250ug

Support: tech@signalwayantibody.com

#AP73396-5 500ug #AP73396-6 1mg

Description

Product Name	Recombinant Human Fibroblast growth factor 1 protein(FGF1) (Active)
Brief Description	Recombinant Protein
Host Species	E.Coli
Uniprot	P05230
GeneID	2246;
Tag Info	NO-tagged
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>

Note: This product is for in vitro research use only