

# Recombinant Rhesus Macaque Interferon gamma protein(IFNG) (Active)



Catalog No: #AP73433

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP73433-1 10ug #AP73433-2 50ug #AP73433-3 100ug #AP73433-4 250ug

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

#AP73433-5 500ug #AP73433-6 1mg

## Description

Product Name	Recombinant Rhesus Macaque Interferon gamma protein(IFNG) (Active)
Brief Description	Recombinant Protein
Host Species	E.Coli
Uniprot	P63310
GeneID	574282;
Tag Info	NO-tagged
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>

Note: This product is for in vitro research use only