

# Recombinant Human Inter-alpha-trypsin inhibitor heavy chain H4(ITIH4),partial



Catalog No: #AP73041

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP73041-1 20ug #AP73041-2 100ug #AP73041-3 1mg

## Description

Product Name	Recombinant Human Inter-alpha-trypsin inhibitor heavy chain H4(ITIH4),partial
Brief Description	Recombinant Protein
Host Species	Yeast
Uniprot	Q14624
GeneID	3700;
Tag Info	N-terminal 6xHis-tagged
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>

Note: This product is for in vitro research use only