

Recombinant Pig Erythropoietin(EPO)

Catalog No: #AP73383



Package Size: #AP73383-1 20ug #AP73383-2 100ug #AP73383-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Pig Erythropoietin(EPO)
Brief Description	Recombinant Protein
Host Species	E.coli
Uniprot	P49157
GeneID	397249;
Tag Info	N-terminal 6xHis-B2M-tagged
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>

Note: This product is for in vitro research use only