

Recombinant Human Inter-alpha-trypsin inhibitor heavy chain H4(ITIH4),partial

Catalog No: #AP71264

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP71264-1 20ug #AP71264-2 100ug #AP71264-3 1mg

Description

Product Name	Recombinant Human Inter-alpha-trypsin inhibitor heavy chain H4(ITIH4),partial
Brief Description	Recombinant Protein
Host Species	E.coli
Uniprot	Q14624
GeneID	3700;
Tag Info	N-terminal 6xHis-SUMO-tagged
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>

Note: This product is for in vitro research use only