

Recombinant Human DNA-directed DNA/RNA polymerase mu(POLM)



Catalog No: #AP77996

Orders: order@signalwayantibody.com

Package Size: #AP77996-1 20ug #AP77996-2 100ug #AP77996-3 1mg

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name | Recombinant Human DNA-directed DNA/RNA polymerase mu(POLM) |
| Brief Description | Recombinant Protein |
| Host Species | Baculovirus |
| Uniprot | Q9NP87 |
| GeneID | 27434; |
| Tag Info | N-terminal 10xHis-tagged and C-terminal Myc-tagged |
| Storage | <p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p> |

Note: This product is for in vitro research use only