

TCTP (phospho-Ser46) rabbit pAb

Catalog No: #13487



Package Size: #13487-1 50ul #13487-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | TCTP (phospho-Ser46) rabbit pAb |
| Host Species | Rabbit |
| Purification | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Applications | WB |
| Species Reactivity | Human,Mouse,Rat |
| Specificity | This antibody detects endogenous levels of Human Mouse Rat TCTP (phospho-Ser46) |
| Immunogen Description | Synthesized phosho peptide around human TCTP (Ser46) |
| Other Names | Translationally-controlled tumor protein (TCTP) (Fortilin) (Histamine-releasing factor) (HRF) (p23) |
| Accession No. | Swiss Prot:P13693GenelD:7178 |
| Uniprot | P13693 |
| GenelD | 7178 |
| SDS-PAGE MW | 24 |
| Concentration | 1 mg/ml |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | -20°C/1 |

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only