

TAZ (phospho-Ser89) rabbit pAb

Catalog No: #13489

Package Size: #13489-1 50ul #13489-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | TAZ (phospho-Ser89) rabbit pAb |
| Host Species | Rabbit |
| Purification | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Applications | WB |
| Species Reactivity | Human,Mouse,Rat |
| Specificity | This antibody detects endogenous levels of Human Mouse Rat TAZ (phospho-Ser89) |
| Immunogen Description | Synthesized phosho peptide around human TAZ (Ser89) |
| Other Names | WW domain-containing transcription regulator protein 1 (Transcriptional coactivator with PDZ-binding motif) |
| Accession No. | Swiss Prot:Q9GZV5GeneID:25937 |
| Uniprot | Q9GZV5 |
| GeneID | 25937 |
| SDS-PAGE MW | 44 |
| Concentration | 1 mg/ml |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | -20°C/1 |

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only