

## TACC3 (phospho-Ser558) rabbit pAb

Catalog No: #13495



Package Size: #13495-1 50ul #13495-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	TACC3 (phospho-Ser558) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of Human TACC3 (phospho-Ser558)
Immunogen Description	Synthesized phospho peptide around human TACC3 (Ser558)
Other Names	Transforming acidic coiled-coil-containing protein 3 (ERIC-1)
Accession No.	Swiss Prot:Q9Y6A5GeneID:10460
Uniprot	Q9Y6A5
GeneID	10460
SDS-PAGE MW	80
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

## Application Details

WB 1:1000-2000

## Background

transforming acidic coiled-coil containing protein 3(TACC3) Homo sapiens This gene encodes a member of the transforming acidic coiled-coil protein family. The encoded protein is a motor spindle protein that may play a role in stabilization of the mitotic spindle. This protein may also play a role in growth a differentiation of certain cancer cells. [provided by RefSeq, Nov 2011],

Note: This product is for in vitro research use only