

SGTA (phospho-Ser305) rabbit pAb

Catalog No: #13530



Package Size: #13530-1 50ul #13530-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SGTA (phospho-Ser305) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of Human SGTA (phospho-Ser305)
Immunogen Description	Synthesized phospho peptide around human SGTA (Ser305)
Other Names	Small glutamine-rich tetratricopeptide repeat-containing protein alpha (Alpha-SGT) (Vpu-binding protein) (UBP)
Accession No.	Swiss Prot:O43765GeneID:6449
Uniprot	O43765
GeneID	6449
SDS-PAGE MW	35
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only