

Semaphorin 4B (phospho-Ser825) rabbit pAb

Catalog No: #13533



Package Size: #13533-1 50ul #13533-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Semaphorin 4B (phospho-Ser825) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	This antibody detects endogenous levels of Human Mouse Rat Semaphorin 4B (phospho-Ser825)
Immunogen Description	Synthesized phosho peptide around human Semaphorin 4B (Ser825)
Accession No.	Swiss Prot:J3KNP4
Uniprot	J3KNP4
SDS-PAGE MW	95
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only