

## SATB1 (phospho-Ser47) rabbit pAb

Catalog No: #13538



Package Size: #13538-1 50ul #13538-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	SATB1 (phospho-Ser47) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	This antibody detects endogenous levels of Human Mouse Rat SATB1 (phospho-Ser47)
Immunogen Description	Synthesized phospho peptide around human SATB1 (Ser47)
Other Names	DNA-binding protein SATB1 (Special AT-rich sequence-binding protein 1)
Accession No.	Swiss Prot:Q01826GeneID:6304
Uniprot	Q01826
GeneID	6304
SDS-PAGE MW	85
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

## Application Details

WB 1:1000-2000

## Background

This gene encodes a matrix protein which binds nuclear matrix and scaffold-associating DNAs through a unique nuclear architecture. The protein recruits chromatin-remodeling factors in order to regulate chromatin structure and gene expression. [provided by RefSeq, Apr 2016].

Note: This product is for in vitro research use only