

RPA32/RPA2 (phospho-Ser8) rabbit pAb

Catalog No: #13550



Package Size: #13550-1 50ul #13550-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RPA32/RPA2 (phospho-Ser8) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of Human RPA32 or RPA2 (phospho-Ser8)
Immunogen Description	Synthesized phospho peptide around human RPA32 and RPA2 (Ser8)
Other Names	Replication protein A 32 kDa subunit (RP-A p32) (Replication factor A protein 2) (RF-A protein 2) (Replication protein A 34 kDa subunit) (RP-A p34)
Accession No.	Swiss Prot:P15927GenelD:6118
Uniprot	P15927
GeneID	6118
SDS-PAGE MW	32
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only