

ROS1 (phospho-Tyr2274) rabbit pAb

Catalog No: #13553

Package Size: #13553-1 50ul #13553-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ROS1 (phospho-Tyr2274) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of Human ROS1 (phospho-Tyr2274)
Immunogen Description	Synthesized phosho peptide around human ROS1 (Tyr2274)
Other Names	Proto-oncogene tyrosine-protein kinase ROS (EC 2.7.10.1) (Proto-oncogene c-Ros) (Proto-oncogene c-Ros-1) (Receptor tyrosine kinase c-ros oncogene 1) (c-Ros receptor tyrosine kinase)
Accession No.	Swiss Prot:P08922GeneID:6098
Uniprot	P08922
GeneID	6098
SDS-PAGE MW	258
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Background

ROS proto-oncogene 1, receptor tyrosine kinase(ROS1) Homo sapiens This proto-oncogene, highly-expressed in a variety of tumor cell lines, belongs to the sevenless subfamily of tyrosine kinase insulin receptor genes. The protein encoded by this gene is a type I integral membrane protein with tyrosine kinase activity. The protein may function as a growth or differentiation factor receptor. [provided by RefSeq, Jul 2008],

Note: This product is for in vitro research use only