

REPS1 (phospho-Ser709) rabbit pAb

Catalog No: #13568



Package Size: #13568-1 50ul #13568-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	REPS1 (phospho-Ser709) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of Human REPS1 (phospho-Ser709)
Immunogen Description	Synthesized phospho peptide around human REPS1 (Ser709)
Other Names	RalBP1-associated Eps domain-containing protein 1 (RalBP1-interacting protein 1)
Accession No.	Swiss Prot:Q96D71GeneID:85021
Uniprot	Q96D71
GeneID	85021
SDS-PAGE MW	80
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Background

RALBP1 associated Eps domain containing 1(REPS1) Homo sapiens This gene encodes a signaling adaptor protein with two EH domains that interacts with proteins that participate in signaling, endocytosis and cytoskeletal changes. The encoded protein has been found in association with intersectin 1 and Src homology 3-domain growth factor receptor-bound 2-like (endophilin) interacting protein 1 when intersectin 1 was isolated from clathrin-coated pits. The encoded protein has also been shown to interact with amphiphysin, a cytoplasmic protein at the surface of synaptic vesicles. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Mar 2014],

Note: This product is for in vitro research use only