

RanBP3 (phospho-Ser58) rabbit pAb

Catalog No: #13577



Package Size: #13577-1 50ul #13577-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RanBP3 (phospho-Ser58) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of Human RanBP3 (phospho-Ser58)
Immunogen Description	Synthesized phospho peptide around human RanBP3 (Ser58)
Other Names	Ran-binding protein 3 (RanBP3)
Accession No.	Swiss Prot:Q9H6Z4GeneID:8498
Uniprot	Q9H6Z4
GeneID	8498
SDS-PAGE MW	60
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Background

This gene encodes a protein with a RanBD1 domain that is found in both the nucleus and cytoplasm. This protein plays a role in nuclear export as part of a heteromeric complex. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008],

Note: This product is for in vitro research use only