

## PPIG (phospho-Ser376) rabbit pAb

Catalog No: #13603

Package Size: #13603-1 50ul #13603-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PPIG (phospho-Ser376) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	This antibody detects endogenous levels of Human Mouse Rat PPIG (phospho-Ser376)
Immunogen Description	Synthesized phosho peptide around human PPIG (Ser376)
Other Names	Peptidyl-prolyl cis-trans isomerase G (PPIase G) (Peptidyl-prolyl isomerase G) (EC 5.2.1.8) (CASP10) (Clk-associating RS-cyclophilin) (CARS-Cyp) (CARS-cyclophilin) (SR-cyclophilin) (SR-cyp) (SRcyp) (Cyclophilin G) (Rotamase G)
Accession No.	Swiss Prot:Q13427GeneID:9360
Uniprot	Q13427
GeneID	9360
SDS-PAGE MW	82
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

## Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only