

PI3 Kinase Class III (phospho-Ser249) rabbit pAb

Catalog No: #13620



Package Size: #13620-1 50ul #13620-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PI3 Kinase Class III (phospho-Ser249) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of Human PI3 Kinase Class III (phospho-Ser249)
Immunogen Description	Synthesized phospho peptide around human PI3 Kinase Class III (Ser249)
Other Names	Phosphatidylinositol 3-kinase catalytic subunit type 3 (PI3-kinase type 3) (PI3K type 3) (PtdIns-3-kinase type 3) (EC 2.7.1.137) (Phosphatidylinositol 3-kinase p100 subunit) (Phosphoinositide-3-kinase class 3) (hVps34)
Accession No.	Swiss Prot:Q8NEB9GeneID:5289
Uniprot	Q8NEB9
GeneID	5289
SDS-PAGE MW	100
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only