

p56Dok-2 (phospho-Tyr351) rabbit pAb

Catalog No: #13645

Package Size: #13645-1 50ul #13645-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	p56Dok-2 (phospho-Tyr351) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of Human p56Dok-2 (phospho-Tyr351)
Immunogen Description	Synthesized phosho peptide around human p56Dok-2 (Tyr351)
Other Names	Docking protein 2 (Downstream of tyrosine kinase 2) (p56(dok-2))
Accession No.	Swiss Prot:O60496GeneID:9046
Uniprot	O60496
GeneID	9046
SDS-PAGE MW	48
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Background

docking protein 2(DOK2) Homo sapiens The protein encoded by this gene is constitutively tyrosine phosphorylated in hematopoietic progenitors isolated from chronic myelogenous leukemia (CML) patients in the chronic phase. It may be a critical substrate for p210(bcr/abl), a chimeric protein whose presence is associated with CML. This encoded protein binds p120 (RasGAP) from CML cells. [provided by RefSeq, Jul 2008],

Note: This product is for in vitro research use only