

Numb (phospho-Ser276) rabbit pAb

Catalog No: #13659

Package Size: #13659-1 50ul #13659-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Numb (phospho-Ser276) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of Human Numb (phospho-Ser276)
Immunogen Description	Synthesized phospho peptide around human Numb (Ser276)
Other Names	Protein numb homolog (h-Numb) (Protein S171)
Accession No.	Swiss Prot:P49757GenelD:8650
Uniprot	P49757
GenelD	8650
SDS-PAGE MW	70
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Background

NUMB, endocytic adaptor protein(NUMB) Homo sapiens The protein encoded by this gene plays a role in the determination of cell fates during development. The encoded protein, whose degradation is induced in a proteasome-dependent manner by MDM2, is a membrane-bound protein that has been shown to associate with EPS15, LNX1, and NOTCH1. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2016],

Note: This product is for in vitro research use only