

## MFF (phospho-Ser146) rabbit pAb

Catalog No: #13709

Package Size: #13709-1 50ul #13709-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	MFF (phospho-Ser146) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	This antibody detects endogenous levels of Human Mouse Rat MFF (phospho-Ser146)
Immunogen Description	Synthesized phosho peptide around human MFF (Ser146)
Other Names	Mitochondrial fission factor
Accession No.	Swiss Prot:Q9GZY8GeneID:56947
Uniprot	Q9GZY8
GeneID	56947
SDS-PAGE MW	37
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

## Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only